hf SURG® Plus

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INHALT

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PLEASE READ THE FOLLOWING INSTRUCTION BEFORE USE.

hf SURG® Plus

EN INSTRUCTIONS FOR USE

CONTENTS

How to use the instructions

Precautions for safety

- Description of product and components
- Appearance and Structure


- Exterior Photos
- Detailed description of the part
- Dimension and Weight
- Composition and Features of the Product
- Summary of EQUIPMENT application specification

- Performance
- Performance

- How to operate or use

HOW TO USE THE INSTRUCTIONS

- Please carefully read the Instructions for the safety before using the product.
- Please see troubleshooting guide in the event of a device failure.
- It's highly recommended to carefully read the Instructions before use to be experienced with operation procedure and precautions without any problems.

- Preparations before use
- How to use
- Practical exercise on a beef model
- Storage and maintenance after use
- How to store and maintain after use
- Troubleshooting
- Operator regulation
- Service
- Disposal of waste products

PRECAUTIONS FOR SAFETY

- ⚠ **Please carefully read the Instructions before use**
 - Do not use near flammable or combustible material to prevent source of potential fire or explosion.
 - Do not disassemble or repair the product without permission. In case of any abnormal operation or malfunction, please stop using the product and contact our office.
 - The replacement parts of the product have to be approved by our company. Warranty of the product will not be accepted in case of using different parts than those supplied by our company.
 - This product may have influence on other electronic equipment and vice versa.
 - Please keep the following Instructions.
- ⚠ **Please do not use the device with the following medical instrument or conditions because it may cause malfunction of the medical instrument. Please do not use together unless otherwise indicated by medical specialist.**

Unless it is recommended by medical specialist

- The patient with electronic implanted medical device.
- The patient with life-support system including respirator.
- The patient with portable electronic medical device including ECG monitor.
- The patient with an electronic implanted device such as a cardiac pacemaker.
- Type of protection against electric shock: Class I equipment
- Degree of protection against electric shock: Type BF applied parts
- Protection against harmful ingress of water (footswitch): IPX8
- Use environment in an Oxygen Rich Environment: Not suitable
- Mode of operation: Operate 10 seconds and wait 30 seconds before repetition
- Operation temperature : 10–35 °C

SYMBOLS USED ON THE PRODUCT

	CE marking with the identification number of the Notified Body (2460 for DNV, Norway) responsible for implementation of the conformity assessment procedures acc. to the MDD (93/42/EEC), Annex V.		
	Symbol for Type BF Applied Parts		Serial no.
	General warning		General cautions
	Date of manufacture		Manufacturer
	Read the instruction		Follow for instruction for use
	Authorized representative in the European community		Crossed-out wheeled bin acc. to WEEE directive (2012/19/EU) for Warning-Disposal of waste products

⚠ Processing instructions according to EN ISO 17664

The product has to be cleaned, disinfected and sterilized according to these processing instructions before its first use as well as after every further use. Disinfection alone is not enough. Thorough cleaning and disinfection are decisive prerequisites for the effective sterilization. The processing should be started as early as possible, 2 hours after use, however, at the latest. During the processing, the product should not be unnecessarily exposed to wetness or humidity. Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the medical practice and/or the hospital.

LIMITATIONS IN THE PROCESSING: The frequent processing has only minor effects on this product. The end of the product's useful life is mainly determined by wear and damage caused by use. In case of doubt, the products should always be sorted out and replaced early. The user alone is responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for the function, performance and safety of the product. **PERSONAL PROTECTIVE EQUIPMENT:** For reasons of occupational safety and to minimize infections / cross infections, suitable personal protective equipment (protective clothes, protective gloves and etc) has to be used during the entire processing.

INSTRUCTIONS: INITIAL TREATMENT AT THE PLACE OF USE: Remove surface contamination using a disposable cloth / paper cloth when you are still at the place of use. Rinse the products with water (at least drinking water quality) 2 hours after the application at the latest. The drying of residues or contamination of any kind on the product has to be avoided. When using the product for the first time at the place of use, do not use any aldehyde-containing or alcohol-containing agents as they may cause protein fixation. **PREPARATION BEFORE CLEANING:** Disassemble the product into its individual components. Visual inspection for damage and wear.

CLEANING: MANUAL The manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified below using a brush below the water surface to achieve both, sufficient cleaning of the products / individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product / individual component; it should, however, not be less than 1 minute. Afterwards, rinsing of the products / individual components with water (at least drinking water quality) for at least 20 seconds. **DESCRIPTION OF THE RECOMMENDED CLEANING PROCESS:** **Pre-cleaning:** Remove surface contamination using a disposable cloth. **Cleaning:** Insert the instrument into the strainer bowl. Hang the strainer bowl into the cleaning bath with cleaning solution, remove contamination by means of cleaning brushes. **Cleaning tools:** Flexbrush brush, REF 605 254, Mirabrush, REF 605 260/61, Apply-Tips, REF 605 510. **Rinsing:** Rinse with water (at least drinking water quality) for 20 sec. **Drying:** Drying at room temperature. The manual cleaning should not exceed a temperature of 45 °C. It has to be ensured that visual contamination has been completely removed. If after the cleaning process, there is still contamination visible, the cleaning has to be repeated.

DISINFECTION: MANUAL For the manual disinfection, please only use approved disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed). Put the products / individual components into corresponding disinfection baths according to the disinfectant manufacturer's specifications. Ensure that the products are sufficiently covered and do not contact each other. The instructions for use of the disinfectant manufacturer must be strictly observed. Particularly the concentrations and contact times have to be observed. Application with the following agents is recommended: - Dürr Dental ID 213 instrument disinfection. **DESCRIPTION OF THE RECOMMENDED DISINFECTION PROCESS:** Completion of proper cleaning preparation and manual cleaning. **Disinfection:** Insert the part of the instrument with patient contact into the strainer bowl. Hang the strainer bowl into the disinfection bath with disinfection solution. **Disinfectant:** ID 213, company Dürr Dental, concentration: 2 %. **Contact time:** 5 min. **Rinsing:** Rinse with water (at least drinking water quality) for 20 sec. **Drying:** Drying at room temperature. Make sure that the disinfectant is compatible with the products and cleaning agents that might be used. The pH of the disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic / halogen hydrocarbons. Please observe thermal restrictions that might apply considering the disinfectant manufacturer's specifications. Afterwards, rinsing of the products / individual components with water (at least drinking water quality) for at least 20 seconds. Before further processing, the product must be free from any residues and dry. **Drying:** If you process the products manually, the products / individual components may be dried by blowing them off with filtered, oilfree compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature.

MAINTENANCE, CONTROL AND CHECK The product does not make special requirements on maintenance. Re-assemble disassembled products / individual components. You must always carry out a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product / individual component is not visibly clean, the entire processing has to be repeated or the product / individual component has to be properly disposed of.

PACKAGING The product is to be packed in suitable and standardized transparent sterile packaging (sterilization bags) and sealed. Observe the instructions of the sterilization bag and sealing machine manufacturers as well as the current, normative requirements. Products / individual parts that have not been sterilized in bags must be used immediately.

STERILIZATION Only tested steam vacuum autoclaves may be used. Make sure that during the sterilization of several different products / individual components, the autoclave is not overfilled and that the products / individual components do not contact each other. The following sterilization cycles may be completed: Steam sterilization, 134 °C, hold time 5 minutes or steam sterilization, 121 °C, hold time 15 minutes To dry the products / individual components, the drying cycle of the autoclave should be set. Use the autoclave manufacturer's instructions for use.

STORAGE To maintain the sterility, the products have to be stored in standardized sterilization bags at a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed once again before their use.

ADDITIONAL INFORMATION Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling. The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded. The manufacturer and the competent authority of your member state shall be immediately notified all serious incidents occurring in connection with the product. For the disposal of the products, there are regionally differing recommendations and regulations. Ask your competent disposal company for the current regulations in your region.

Cleaning and Sterilizations

HANDPIECE, RETURN ELECTRODE, FOOTSWITCH, POWER CORD: Always clean after each use patient. Remove the electrode from the handpiece to cleaning. The outer surface of the handpiece, return electrode, footswitch and power cord should be cleaned with antiseptic soap or solution; wipe with disinfectant and/or water. **ELECTRODES:** Always clean and sterilize after each use to prevent spreading germs and diseases to other patients. Saliva, blood, and other debris may be left on the electrode. To sterilize electrode, wipe the electrode thoroughly using a mild detergent. Rinse thoroughly and dry. Steam sterilizing with 134 °C, 15 minutes, 2.1kgf/cm²


Cleaning and Sterilizations

DEVICE. Since the device does not have direct contact with the patient, it can be cleaned by carefully wiping with alcohol. Avoid using disinfectants not specified for use with metal.

The device is used by medical specialist and doctor. In case of the followings, DNHmedtech will not responsible for any loss which may result in the user's wrong use.

- Do not tamper with device or any of its components.
- For safety operation, confirm that the device is plugged into a three-wire grounded power source. Do not use plug adaptors not ground.
- Do not use on patient with active implants such as pacemakers, cochlear or nerve stimulators.
- Metal can conduct radio frequency energy just same as electricity. Remove partial dentures if it is removable and check the area around fixed partial dentures. Metal restorations can conduct energy and cause burns and discomfort to the patient.
- Be careful that metal restorations, bone or teeth do not come into contact with active electrodes. While transient contact is to be detrimental, extended contact may result in damage. All instruments used during device should be made of a non-conducting plastic material. Cotton rolls or gauze pads used in the patient's mouth should be kept from moist during operation. Skin-to-skin contact (for example between the arms and body of the patient) should be avoided by insertion of dry gauze for an example.
- When high frequency surgical and physiological monitoring equipment are used simultaneously on the same patient, all monitoring equipments of electrodes should be placed as far as possible from the surgical electrode.
- Position the cables to the surgical electrodes in such a way that contact with the patient or so that other leads can be avoided. Store electrodes out of reach the patient.
- Return electrode (patient plate) must be used during all procedures. The entire area of return electrode should be placed firmly, nonconductive contact (clothes on) with the patient, preferably against the patient's upper back, contacting the maximum possible area.
- Do not place the return electrode (patient plate) near the heart.
- The use of flammable anesthetics or oxidizing gases such as nitrous oxide and oxygen should be avoided. Non-flammable agents should be used for cleaning or disinfecting. Solvents of adhesives is allowed to evaporate before the application. Some material may be ignited by sparks produced in normal use of the equipment (for example, cotton wool and gauze when saturable with oxygen). Gases within the device may be ignited. The device should be used after completely removing anesthetics or disinfections.
- The lowest effective output setting could be achieved the best results.
- The interference produced by the operation of the high frequency may adversely influence on the operation of other electrical equipment. In case of interference, de-energize or increase the distance between equipments. Connecting with a different power circuit may also reduce interference.
- Regularly inspect the accessories, particularly the electrode sheaths and cables, to avoid possible insulation damage.
- The electrode is composed with the tip (wire, loop, ball etc) and insulation part. The incision and coagulation have to be operated with tips of electrode. The electrode insulation part have a limited life expectancy and should be replaced regularly or exchanged according to the deterioration condition.
- Be sure to use appropriate accessories supplied from DNHmedtech. Unauthorized cables and accessories may negative effect on EMC performance and safe operation.
- Failure of the device or excessive electrical interference could result in an unintended power increase, decrease or activation. In case of electrical interference, de-energize or increase the distance between equipments causing interference. Connecting to a different power circuit may also reduce interference.
- Recommended not to use the type of needle-shaped monitoring electrodes.

1. DESCRIPTION OF PRODUCT AND COMPONENTS

No.	Name	Doc No.	QTY	Description
1	Device	ES-SET-50A	1	Electrosurgical Unit Monopolar
2	Handpiece(assembly)	ES-050-300	1	HF energy of the electrosurgery travels from the monopolar handpiece and electrode tip to return electrode.
3	Handpiece table	ES-050-304	1	Handpiece holder for resting
4	Return electrode (patient plate / assembly)	ES-050-301	1	HF energy of electrosurgery returns to return electrode.
5	Electrode tip	a D15-45	1	Coagulation: Ball electrode Ø1.6 mm shaft 45° angled
		b D01-45	1	Knife: Needle electrode 0.2mm straight shaft 45° angled, L 7.0 mm
		c D03-45	1	Scoop: Wire loop electrode elongated 6 x 1.7 mm shaft 45° angled, wire 0.2 mm
		d D12-90	1	Coagulation: Coagulation and fulguration electrode 0.5 mm 90° angled, L 10 mm
		e D05-45	1	Scoop: Wire loop electrode elongated 6 x 1.7 mm, 45° angled, wire 0.2 mm
		f D06-00	1	Scoop: Wire loop electrode Ø 4 mm straight, wire 0.2 mm
6	Footswitch(assembly)	ES-050-303	1	Determine the performance level
7	Power cord	ES-050-306	1	The AC main power supply
8	User manual	ES-050-307	1	Operation and instruction manual

2. APPEARANCE AND STRUCTURE

1. EXTERIOR PHOTOS

1) FRONT VIEW

2) REAR VIEW


3) SIDE VIEW (RIGHT)

(LEFT)

2. DETAILED DESCRIPTION OF THE PART

THE FRONT

- a) Power switch: Main power on or off
- b) Power lamp: When power switch (a) is on, light turns green.
- c) Function (CUT) switch: Initially when power switch is on, start function "CUT", by pressing button.
- d) Function (CUT) lamp: If "CUT" is selected, the yellow light is on.
- e) Function (CUT/COAG): Selected function by pressing button
- f) Function (CUT/COAG) lamp: If "CUT/COAG" is selected, the yellow light is on.
- g) Function (COAG) switch: Selected function by pressing button
- h) Function (COAG) lamp: If "COAG" is selected, the blue light is on.
- i) Power rotary switch: Select output power level (1~10).

- j) Power level lamp: Power levels are displayed cumulatively of green lights by power rotary switch
- k) RF/Output lamp: The orange light is on when the output is operated by footswitch.
- l) Handpiece socket
- m) Return electrode (patient plate) socket
- n) Footswitch socket

THE REAR

- a) AC power inlet: Connect AC power cord
- b) Fuse cartridge: Main fuses (250V, F2.5AHx2)
- c) Volume: Tuning of the tone, when operated.
- d) Label

3. DIMENSION AND WEIGHT

1) Dimension of device: approx. (W) 265 x 180 x 89 mm 2) Weight of device (only): approx. 3.5 kg 3) Dimension and weight of component

Name	Appearance	Dimension	Weight
Handpiece		Handle length: 15 cm Wire: 200 cm	50 g
Return electrode (Patient plate)		Width: 9cm Height: 9cm Wire: 200cm	50 g
Electrode set		Length: 45mm	6 g
Footswitch		Width: 7 cm Height: 3 cm Length: 10 cm Wire: 200 cm	290 g
Power Cord		150 cm	210 g

4. COMPOSITION AND FEATURES OF THE PRODUCT

1) Product composition: The device is used to generate high frequency for cutting and coagulation of tissue. It is composed of active electrode (Handpiece), return electrode (Patient plate), electrode, foot switch and the AC power cord supplying the main power.

2) Product characteristic:

- a) Characteristic : The device using generated high frequency energy is able to cut and coagulate tissue with focused heat energy at the active electrode-tip. It has three output modes: 'CUT', 'CUT&COAG' and 'COAG'. The intensity of these modes and output levels can be adjusted by the operator in accordance with purpose of treatment.
- b) Operation principle: The high frequency energy is generated at the main device, it is returned to electrode through the patient's body placed on the patient's back. At this time high frequency energy is created at the electrode and hyper-thermal energy is placed at the treatment area. Depending on the intended use, the output energy of device may be either partially or fully rectified. A fully rectified output is appropriate for cutting, while coagulation is achieved with a partially rectified output.
- c) The electrical specification
 - The rated voltage: AC 220-230V
 - The rated frequency: 50/60Hz
 - Power consumption : Actto-50A (220VA)
- d) The protection type and degree for the electric shock: the first class equipment, BF type equipment.
- e) The safety device
 - The power transformer is blocked by the fuses inside equipment when an over current is detected.
 - The output and power levels are controlled by output level switch and foot switch in accordance with treatment purpose. And operator or staff

can check all device conditions such as function lamps ("CUT": yellow, "CUT/COAG": yellow, "COAG": blue), output power level lamp (green) and operation lamp ("RF/OUTPUT": orange)

- The device emits a tone while the output is energized.

5. SUMMARY OF EQUIPMENT APPLICATION SPECIFICATION

- 1) Medical purpose. Electrosurgical unit, monopolar type Actto-50 is a high frequency electrosurgical system, intended to cut and coagulation tissue by the radio frequency energy, focusing the heat energy at the small and active electrode.
- 2) Appliessd part of body or type of tissue interacted with
 - a) Treatment site: teeth areas
 - b) Condition : intact teeth areas within oral orific
- 3) Intended USER
 - a) Language understanding: - language as specified in the operation manual
 - b) Experience: - no special experience needed - no maximum
 - c) Permissible impairments:
 - mild reading vision impairment or vision corrected tolog MAR 0.2 (6/10 or 20/32);
 - one arm/hand system capable of guiding and holding device;
 - average degree of aging-related short term memory impairment;
 - impaired by 40 % resulting in 60 % of normal hearing at 500 Hz to 2 kHz
- 4) Environment
 - a) General: - hospital, intended for professional use only;
 - indoor use only;
 - b) Conditions of visibility:
 - Ambient luminance 100-500 lux
 - Viewing distance 20 cm to 1 metre
 - Viewing angle: normal to the scale $\pm 20^\circ$
 - c) Physical
 - temperature range: 10 °C to 35 °C
 - relative humidity range: 30 % to 75 %
 - Atmospheric pressure range: 70 to 106 kPa
 - background sound pressure level: <70 dB(A) in the range of 100 Hz - 8 kHz
 - d) Frequency of use: - unlimited; but 10s On and 30s Off
 - e) Mobility: - transportable; but placed on the table during use

3. PERFORMANCE

1. PERFORMANCE

- 1) Protection type and level for electric shock : Class1 Equipment, Type BF Applied Part
- 2) The power supply device
 - Input voltage and power consumption AC 220-230V, 50/60 Hz, Actto-50A(220VA)
 - Power transformer output voltage : Actto-50A (83V, 7V)
- 3) The output
 - a) Normal operating frequency : Actto-50A (1.2 MHz $\pm 10\%$)
 - b) Coagulation modulation wave shape : Square wave
 - c) Maximum output power (Output impedance : Actto-50A (400 Ω))
 - CUT : 50W $\pm 20\%$ - CUT/COAG : 45W $\pm 20\%$ - COAG : 30W $\pm 20\%$
 - d) Maximum output voltage (Output impedance : Actto-50A (400 Ω))
 - CUT : 280 Vpeak - CUT/COAG : 270 Vpeak - COAG : 260 Vpeak
 - e) Rated accessory voltage - Handpiece : 600 Vpeak - Return electrode : 600 Vpeak

4. HOW TO OPERATE OR USE


1. PREPARATIONS BEFORE USE

- 1) Read the instruction manual to get used to operation and precautions before use.
- 2) Check all components.
- 3) Before connecting AC power supply, check that the power switch is in the "Off" position.
- 4) Place the device in a location where it will not be potentially explosive atmospheres because electrical spark may occur within the device.
- 5) Do not leave the device in wet location.
- 6) Check if the components (Handpiece, Return electrode, Footswitch) are connected correctly.
- 7) When necessary, clean and sterilize electrode sheath by a gravity steam sterilization at 134 °C for 15 minutes.
- 8) Combine the appropriate electrode into the handpiece holder, fix it firmly to the handpiece by turning it clockwise. Check the tightness of the electrode by gently pulling it. It has to be impossible to turn the electrode or to drag it out of the shaft. Make sure that the non-coated part of the electrode (metal part) is inserted completely in the handpiece. No metal parts of the shaft should be visible outside of the head of the handpiece. Otherwise, you pose a risk of injury to you and your patients. notice: Only the approved original electrodes must be used. These electrodes are adapted to the handpieces and devices.

2. HOW TO USE

- 1) Connect AC Power cord to power inlet at the back of device. And connect the plug to a grounded power outlet. (The appliance coupler shall be placed to be easily accessible.)
- 2) Turn the power switch on ("I"), and check power lamp is on(green).
- 3) Select the functions ("CUT";"CUT/COAG";"COAG") by pressing the button, and check display lamp.
- 4) Fasten the return electrode on patient's back or insert the patient plate to patient's waist belt. Notice: When operating the electrosurgery, it must always be worked with the connected return electrode. This ensures an optimal power.
- 5) Before applying the electrode to patient tissue, check the normal operation by the OUTPUT lamp (orange) and the acoustic tone by pressing the footswitch. Adjust volume of the alarm.
- 6) Set output power level by power rotary switch and check of green lights.
- 7) The lowest output set is effective for the best results. If the electrode is effected on electrical resistance, the set level is too low. Increase the level until there is no resistance for cutting, sparking and discoloration along the incision. This point is just the lowest effective value.
- 8) Operating time is repeated at 30 seconds intervals over the following 10 seconds.
- 9) Position cables to surgical electrodes in such a way that they are contacted with the patient or other leads are not contacted with the patient. Store active electrodes out of reach the patient.

Note: After reading instruction, practicing on raw and lean beef can help the practitioner acquire the necessary skills to achieve superior clinical results.

3. PRACTICAL EXERCISE ON A BEEF MODEL

Prepare the device for operation and follow the steps mentioned below.

- 1) Select a piece of fresh, lean beef. Veal is not suitable, because it does not change color when cut with an electrode. Because of its cell structure pork is also not suited. Wait until the beef has obtained room temperature. Note: Make sure that the meat is placed on the plugged-in neutral electrode. If this is not the case the waves cannot derivate and, consequently cannot proceed with the exercise.
- 2) Insert the electrode of your choice (Knife tip, loop etc) into the hand piece.
- 3) Turn the intensity control (power) dial to 10
- 4) Push the function button to "CUT"
- 5) Activate the foot switch
- 6) Make several incisions of different lengths and depths with even, brushing movements. Then take the energy from the electrode and look at the result. You will notice that the intensity adjustment was too high, which caused sparks and remarkable discoloration along the cutting line.
- 7) Reduce the intensity setting to 1 (or low power). You will notice that the electrode will either cut only if it is dragged and pulled through the meat, or does not cut at all. If a cut is possible at all, bits and pieces of meat will get caught on the electrode.
- 8) Repeat the above described process with slightly increased settings, until you reach the point at which discolorations and visible discharges of sparks no longer appear. The tip of the electrode should not encounter resistance. The cut should be precisely even, without the occurrence of a discharge of sparks and without the necessity to drag the electrode. Continue with your endeavors by executing slow, intermediate and faster cuts for each particular setting in order to achieve the necessary expertise and confidence you will need for the surgery of your patients.
- 9) Turn the function switch to cut with simultaneous coagulation and repeat the exercises. You will see that cuts with the slightly modified wave require a higher setting than when the fully filtered wave is used. This is normal and should be taken into consideration for the later work on the patient.

APPLICATION EXAMPLES: HF

Notice: For reference only, the grade should be adjusted according to the user's level and procedure method.

Function	Grade	Power(Watt)	Indications & remarks
CUT	6	27	• sulcus dilatation • internal gingivectomy • tumor resection • gingivectomy • open curettage • lap preparation • excision
CUT	5	21	• vestibulum plastic
CUT/COAG	6	22	• gingivoplastic • exposure of teeth, stubs, approximal steps or crown edges • removal of hyperplasia for ablation of tissue if simultaneous coagulation is requested with the cut (notice: 10% loss of tissue about 24 h post operatively due to extended lateral heat)
CUT/COAG	5	17	Use only if distance to bone or periost is sufficient!

5. STORAGE AND MAINTENANCE AFTER USE

1. HOW TO STORE AND MAINTAIN THE PRODUCT AFTER USE

- 1) After using the device, the power switch should be disconnected.
- 2) When storing the accessories (Power cord, Handpiece, Return electrode, Footswitch, Electrode) separate them from the device, do not force it and check the disconnected cables are not twisted. Especially, keep return electrode (patient plate) packaged with wire and plate at all times.
- 3) If necessary, clean and sterilize electrode sheath by a gravity steam sterilizing at 134 °C for 15 minutes.
- 4) The device should be kept clean for next use. Wipe the exterior surface with a soft cloth and some neutral detergent.
- 5) The product should be stored at the place:
 - In a dry condition
 - Away from direct sunshine, dust, salt or sulfur which may have negative effect on surrounding environment.

	Atmospheric pressure	Temperature	Humidity
condition	70~106 kPa	+10~+35 °C	30~75 %

- Storage and transportation -20 to +60 °C, 15 to 93 % r.h., 50 to 106 kPa

2. TROUBLESHOOTING

- 1) In the event of power failure, check the power cord, connection (plug and socket) and fuses.
- 2) If the temperature is not controlled, check the instruction manual.
- 3) If there is no output power under normal conditions (RF lamp on: orange, BEEP on), check the connection with handpiece and return electrode (patient plate).

3. OPERATOR REGULATION

The device unit is classified as medical device unit class IIb (Europe).

Technical controls: The user is committed to perform on a regular basis technical controls after the following specifications.

Period: every 24 months, starting with date of delivery and after each repair.

Covering: Visual check of the unit and accessories - Check according to IEC 62353

- Protective earth resistance

- Alternate leakage current

- Alternate patient leakage current

Function check:

- Main switch

- Switch CUT, CUT/COG, COAG with LEDs

- Uniformity of performance throughout the range of adjustment

Measurement of radio frequency (HF) output power at a load of 400Ω:

- Output: CUT (50 W ± 20 %), CUT/COAG (45 W ± 20 %), COAG (30 W ± 20 %)

All results of measurements must be documented according to DIN EN ISO 62353 concerning the first measured values. If defects occur during the controls, the user is responsible to initiate repair.

4. SERVICE

If requested, circuit diagrams, components part lists, descriptions, calibration instructions or other information is provided for service activity of the equipment from DNHmedtech Co.

5. DISPOSAL OF WASTE PRODUCTS


When the crossed-out wheeled bin symbol is attached to a product, it means the product is followed by the

European Directive 2012/19/EU. All electrical and electronic products should be disposed of separately to

designated collection facilities of the municipal waste appointed by the local government or authorities.

The correct disposal of old appliance will help to prevent potential negative consequences for the environment and human health. For more detailed information about disposal of waste products, please contact the local authorities or waste disposal service.

Note: Please contact the supplier or the head office for inquiries.

Hager & Werken GmbH & Co. KG

Ackerstraße 1, 47269 Duisburg, Germany

Tel.: +49 (203) 99 26 9-0, Fax: +49 (203) 29 92 83

www.hagerwerken.de

CALMED INVEST Kft.(HU11678179)

1191 Budapest, Attila ucta 34., Budapest, Hungary

Phone)+36 1 291 5066 Fax)+36 1 291 0366

DNHmedtech Co.

#302 B-dong, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si,
Gyeonggi-do, Korea

www.DNHmedtech.com (e-Mail: DNHmedtech@naver.com)

Phone)+82 70 4242 1011 Fax)+82 50 4218 1164

REF 452 500

WEEE NR. DE21760541

CE 2460

2024/42/2013 Rev.2 (2024-06)