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Intended use
Dental splinting material long-term

Extended purpose
Tough-elastic, highly transparent thermoforming foil for the production of occlusal splints such as orthodontic devices or aligner and retainer splints.

Indications for use
Due to its properties, the material is suitable for the production of bite and bruxism splints, as well as for the treatment of myoarthropathies (muscle/jaw joint disorders) and resulting parafuncions such as clenching, grinding etc.. In orthodontics, Biolon can also be used to treat malocclusions.

Contraindication
In case of allergies to any of the ingredients do not use Biolon. Not desired side effects of this medical device are not to be expected in case of proper processing and application. However, immune reactions (e. g. allergies) or irritations principally cannot be excluded. Should you get to know non-desired side effects – also in doubtful cases – we appreciate in any case to be informed with a detailed description of the attendant circumstances and symptoms. We will follow each hint.

Target patient group
Patients who need dental care.

Intended users
Trained staff in dental field

Processing
All Biolon materials are delivered with an insulating foil. This thin foil (0.07 mm) is thermoformed together with the Biolon on the model side and guarantees a decisively improved transparency of the splint. Thus an additional insulation of the model surface is no longer necessary.

If processed in a Vacfomat, the insulation foil must be removed before thermoforming.

The adjustment of the heating and cooling times can either be taken from the survey of thermoforming foils, or can be programmed directly into the unit with the help of the barcode on the Biolon packaging.

If a Drufoformat scan or a Drufosmart scan is used, please take note of the special information included in this package. We recommend Cenit L (REF D3815) for the necessary blocking out works that are done on the model. The splints produced from Biolon can be built up with dental acrylic resins.

Cleaning
For cleaning, the use of MaxiFresh™ cleaning tablets is recommended.
Please be aware of the MaxiFresh™ instruction for use.

Note
While working out on the material, it is recommended to use an exhaust system and protective glass or goggles as well as a face mask.

Serious incidents
All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

Storage conditions
At room temperature, i. e. 18–28 °C.
Protect against humidity, always close foil bag when packs are in use.

Lot number / Expiry date
The lot number and the expiry date are indicates on each Biolon packaging. In case of claims, please always indicate the lot number of the product. Do not use the product after the expiry date.

Disposal
Disposal of contents / the packaging in accordance with local / regional / national / international regulations.

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This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product

according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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