

Technical Data Sheet

Product description: Brushable separating agent for insulating silicones and polyether.

Manufacturer: Drewe Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical designation
1	distilled water
2	polyvinyl alcohol
3	coating additives

Properties:

Color: achromatic
Form: liquid
Odor: characteristic
Boiling point: 100°C
Viscosity: $\leq 1 \text{ Pa}\cdot\text{s}$

Drewe Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

dentamid@drewe.de

www.drewe.com