

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 1 of 8

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

VITA ADIVA CERA ETCH

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

###### Use of the substance/mixture

Use as laboratory reagent

##### 1.3. Details of the supplier of the safety data sheet

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
e-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
e-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	

##### 1.4. Emergency telephone number:

+49-(0)761-19240

#### SECTION 2: Hazards identification

##### 2.1. Classification of the substance or mixture

###### UN GHS (ST/SG/AC.10/11/Rev.8)

Hazard categories:

Acute toxicity: Acute Tox. 2

Acute toxicity: Acute Tox. 3

Acute toxicity: Acute Tox. 4

Skin corrosion/irritation: Skin Corr. 1

Serious eye damage/eye irritation: Eye Dam. 1

Hazard Statements:

Fatal in contact with skin.

Toxic if swallowed.

Harmful if inhaled.

Causes severe skin burns and eye damage.

Causes serious eye damage.

##### 2.2. Label elements

###### UN GHS (ST/SG/AC.10/11/Rev.8)

###### Hazard components for labelling

hydrofluorid acid

Sulphuric acid ... %

Signal word: Danger

Pictograms:



## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 2 of 8

#### Hazard statements

H301	Toxic if swallowed.
H310	Fatal in contact with skin.
H314	Causes severe skin burns and eye damage.
H332	Harmful if inhaled.

#### Precautionary statements

P260	Do not breathe dust/fume/gas/mist/vapours/spray.
P405	Store locked up.
P308+P311	IF exposed or concerned: Call a POISON CENTER/doctor.
P280	Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

#### 2.3. Other hazards

No information available.

### SECTION 3: Composition/information on ingredients

#### 3.2. Mixtures

##### Hazardous components

CAS No	Chemical name	Quantity
	Classification (UN GHS (ST/SG/AC.10/11/Rev.8))	
7664-93-9	Sulphuric acid ... %	5 - < 10 %
	Skin Corr. 1A; H314	
7664-39-3	hydrofluorid acid	1 - < 5 %
	Acute Tox. 1, Acute Tox. 2, Acute Tox. 2, Skin Corr. 1A; H310 H330 H300 H314	
64-17-5	ethyl alcohol	1 - < 5 %
	Flam. Liq. 2, Eye Irrit. 2; H225 H319	

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### General information

First aider: Pay attention to self-protection! Remove affected person from the danger area and lay down.

##### After inhalation

Provide fresh air. Medical treatment necessary. IF INHALED: Call a POISON CENTER/doctor// if you feel unwell.

##### After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Call a physician immediately. IF ON SKIN (or hair): Ca-Gluconate solution  
Immediately call a POISON CENTER/doctor/.

##### After contact with eyes

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.Immediately call a POISON CENTER/doctor/.

##### After ingestion

Rinse mouth immediately and drink plenty of water. Do NOT induce vomiting. Adverse human health effects and symptoms: Gastric perforation. Call a physician immediately. If swallowed, immediately drink: Ca-Gluconate solution

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 3 of 8

#### **4.2. Most important symptoms and effects, both acute and delayed**

No information available.

#### **4.3. Indication of any immediate medical attention and special treatment needed**

Treat symptomatically.

### SECTION 5: Firefighting measures

#### **5.1. Extinguishing media**

##### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

#### **5.2. Special hazards arising from the substance or mixture**

Non-flammable.

#### **5.3. Advice for firefighters**

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

#### **Additional information**

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

### SECTION 6: Accidental release measures

#### **6.1. Personal precautions, protective equipment and emergency procedures**

##### **General advice**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

#### **6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

#### **6.3. Methods and material for containment and cleaning up**

##### **Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### **6.4. Reference to other sections**

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

### SECTION 7: Handling and storage

#### **7.1. Precautions for safe handling**

##### **Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

##### **Advice on general occupational hygiene**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

#### **7.2. Conditions for safe storage, including any incompatibilities**

##### **Requirements for storage rooms and vessels**

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Provide adequate ventilation as well as local exhaust at critical locations.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 4 of 8

#### Hints on joint storage

No special measures are necessary.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Exposure limits (EH40)

CAS No	Substance	ppm	mg/m <sup>3</sup>	fibres/ml	Category	Origin
64-17-5	Ethanol	1000	1920		TWA (8 h)	WEL
7664-39-3	Hydrogen fluoride (as F)	1.8	1.5		TWA (8 h)	WEL
		3	2.5		STEL (15 min)	WEL
7664-93-9	Sulphuric acid (mist)	-	0.05		TWA (8 h)	WEL

### 8.2. Exposure controls



#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

Recommended glove articles KCL Dermatril P Wearing time with occasional contact (splashes): 480 min When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

##### Skin protection

Wear suitable protective clothing.

##### Respiratory protection

If handled uncovered, arrangements with local exhaust ventilation have to be used.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state:

Colour: light red

Odour: characteristic

#### Changes in the physical state

Melting point/freezing point: not determined

Boiling point or initial boiling point and boiling range: 100 °C

#### Flammability

Solid/liquid: not applicable

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 5 of 8

Gas:	not applicable
Lower explosion limits:	not determined
Upper explosion limits:	not determined
<b>Self-ignition temperature</b>	
Solid:	not applicable
Gas:	not applicable
Decomposition temperature:	not determined
pH-Value:	2,0
<b>Solubility in other solvents</b>	
not determined	
Partition coefficient n-octanol/water:	not determined
Vapour pressure: (at 50 °C)	<=1100 hPa
Density:	1,06000 g/cm <sup>3</sup>
Relative vapour density:	not determined

#### 9.2. Other information

##### Information with regard to physical hazard classes

Oxidizing properties

Not oxidising.

##### Other safety characteristics

Solid content: 1,5 %

Evaporation rate: not determined

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

#### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

#### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

#### 10.4. Conditions to avoid

none

#### 10.5. Incompatible materials

No information available.

#### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

### SECTION 11: Toxicological information

#### 11.1. Information on hazard classes as defined in GB CLP Regulation

##### Acute toxicity

Fatal in contact with skin.

Toxic if swallowed.

Harmful if inhaled.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 6 of 8

#### ATEmix calculated

ATE (oral) 100,4 mg/kg; ATE (dermal) 100,4 mg/kg; ATE (inhalation vapour) 10,04 mg/l; ATE (inhalation dust/mist) 1,004 mg/l

CAS No	Chemical name				
	Exposure route	Dose	Species	Source	Method
7664-39-3	hydrofluorid acid				
	oral	ATE 5 mg/kg			
	dermal	ATE 5 mg/kg			
	inhalation vapour	ATE 0.5 mg/l			
	inhalation dust/mist	ATE 0.05 mg/l			
64-17-5	ethyl alcohol				
	dermal	LD50 7060 mg/kg			

#### Irritation and corrosivity

Causes severe skin burns and eye damage.

Causes serious eye damage.

#### Sensitising effects

Based on available data, the classification criteria are not met.

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP]. Special hazards arising from the substance or mixture!

## SECTION 12: Ecological information

### 12.1. Toxicity

The product is not: Ecotoxic.

### 12.2. Persistence and degradability

The product has not been tested.

### 12.3. Bioaccumulative potential

The product has not been tested.

### 12.4. Mobility in soil

The product has not been tested.

### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

### 12.7. Other adverse effects

No information available.

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 7 of 8

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Disposal recommendations

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

#### Contaminated packaging

This material and its container must be disposed of as hazardous waste. Handle contaminated packages in the same way as the substance itself.

## SECTION 14: Transport information

### Marine transport (IMDG)

<b>14.1. UN number or ID number:</b>	UN 2922
<b>14.2. UN proper shipping name:</b>	CORROSIVE LIQUID, TOXIC, N.O.S. (Wasser)
<b>14.3. Transport hazard class(es):</b>	8
<b>14.4. Packing group:</b>	II
Hazard label:	8+6.1



Special Provisions:	274
Limited quantity:	1 L
Excepted quantity:	E2
EmS:	F-A, S-B

### Air transport (ICAO-TI/IATA-DGR)

<b>14.1. UN number or ID number:</b>	UN 2922
<b>14.2. UN proper shipping name:</b>	CORROSIVE LIQUID, TOXIC, N.O.S. (Wasser)
<b>14.3. Transport hazard class(es):</b>	8
<b>14.4. Packing group:</b>	II
Hazard label:	8+6.1



Special Provisions:	A3 A803
Limited quantity Passenger:	0.5 L
Passenger LQ:	Y840
Excepted quantity:	E2
IATA-packing instructions - Passenger:	851
IATA-max. quantity - Passenger:	1 L
IATA-packing instructions - Cargo:	855
IATA-max. quantity - Cargo:	30 L

### 14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

### 14.6. Special precautions for user

Warning: Toxic. strongly corrosive.

### 14.7. Maritime transport in bulk according to IMO instruments

not applicable

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA CERA ETCH

Revision date: 12.02.2020

Product code: 307

Page 8 of 8

#### SECTION 15: Regulatory information

##### National regulatory information

Employment restrictions:	Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers.
Skin resorption/Sensitization:	Permeates easily through outer skin and causes poisoning.

#### SECTION 16: Other information

##### Changes

This data sheet contains changes from the previous version in section(s): 1,2,4,6,7,8,9,10,11,14,15.

##### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%

##### Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
Printing date: 27.02.2017

Version: 1

### 1. Identification of the substance/mixture and of the company/undertaking

#### Product Identifier

Trade name: VITA ADIVA C-PRIME

Use of the preparation: Dental product

Details of the supplier of the safety data sheet:

#### Manufacturer/supplier:

Harvard Dental International GmbH,  
Margaretenstraße 2-4,  
15366 Hoppegarten, Germany  
info@harvard-dental-international.de

Emergency telephone number:

+ 49 (0) 30 99 28 978-0 (Quality Management)

### 2. Hazards identification

Classification of the substance or mixture:

Classification according to Regulation (EC) No 1272/2008:

Flam. Liq. 2, H225 Highly flammable liquid and vapour

Information concerning particular hazards for human and environment:

The product has to be labeled due to calculation procedure for the "General Classification guideline for preparations of the EU" in the latest version.

Label elements

Labeling according to Regulation (EC) 1272/2008:

Hazard pictograms



Signal word

Danger

Hazard statements

H225 Highly flammable liquid and vapour

Precautionary statements

P210 Keep away from heat/sparks/open flames/hot surfaces. – No smoking.

Labeling according to EU guidelines:

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

Other hazards

Results of PBT and vPvB assessment

PBT: not applicable

vPvB: not applicable

## Safety Data Sheet

### according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

<b>3. <u>Composition/Information on Ingredients</u></b>
<p><b>Chemical characterization:</b>  <b>Description:</b>          Solution of methacrylsilanes in ethanol  <b>Hazardous components:</b>          64-17-5 (CAS-No.)                  Ethanol                  GHS02, H225</p>
<b>4. <u>First-aid measures</u></b>
<p><b>General information:</b> Seek medical advice, if symptoms occur, which can be caused by the product  <b>after skin contact:</b> wash with plenty of water and soap  <b>after inhalation:</b> fresh air  <b>after eye contact:</b> Rinse opened eye for several minutes under running water and consult a doctor  <b>after swallowing:</b> Consult a physician.</p>
<b>5. <u>Fire fighting measures</u></b>
<p><b>Extinguishing media:</b>          CO<sub>2</sub>, powder, water spray, foam.          Use extinguishing method and media suitable to surrounding fire.  <b>Special hazards arising from the substance or mixture</b> No further relevant information available  <b>Protective equipment</b> Wear self-contained respiratory protective equipment</p>
<b>6. <u>Accidental release measures</u></b>
<p><b>Personal precautions, protective equipment and emergency procedures:</b>          Do not inhale vapours/aerosols. Avoid substance contact. Ensure supply of fresh air in enclosed rooms  <b>Environmental precautions:</b>          Do not allow to enter drains/surface water/ ground-water.  <b>Methods for cleaning up:</b>          Absorb with liquid binding material. Clean up affected area.  <b>Reference to other sections:</b>          See Section 7 for information on safe handling          See Section 8 for information on personal protection equipment          See Section 13 for disposal information</p>
<b>7. <u>Handling and storage</u></b>
<p><b>Handling:</b>          This product should only be supplied or handled to dentists or dental technicians or upon their instructions.  <b>Precautions for safe handling:</b> Observe normal care for working with chemicals  <b>Information about fire and explosion protection:</b> Keep away from heat and fire. Take measure to prevent electrostatic charge  <b>Storage:</b> Store at a dry place. Store tightly closed. Storage temperature &lt;25°C.  <b>Requirements at storerooms and containers:</b> Store only in the original containers  <b>Further information on storage conditions:</b>          Store away from food stuffs.  <b>Specific end uses:</b> No further relevant information available</p>

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

### 8. Exposure controls/personal protection

**Components with limits of values to be supervised at the workplace:**

CAS-No.	name	art	value	unit
64-17-5	Ethanol	TRGS 900	500	ml/m <sup>3</sup> (ppm)
			or 960	mg/m <sup>3</sup>

**Additional information:** The lists valid during the making were used as basis

**Personal protective equipment:**

**General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

Keep away from foodstuffs and beverages. Avoid contact with eyes or skin

**Hand protection**

Protective gloves.

Due to missing tests no recommendation to the glove material can be given for the product.

**Eye protection:** safety glasses

**Respiratory protection:** Required, when aerosols are generated

### 9. Physical and chemical properties

**Information on basic physical and chemical properties**

**Appearance**

<b>Form:</b>	liquid			
<b>Color:</b>	colorless			
<b>Odor:</b>	characteristic			
	<table border="1"> <thead> <tr> <th>Value</th> <th>Unit</th> <th>Method</th> </tr> </thead> </table>	Value	Unit	Method
Value	Unit	Method		

**Safety relevant basic data**

<b>Melting point/melting range:</b>	-114	°C
<b>Boiling point/boiling range:</b>	78	°C.
<b>Flash point:</b>	12	°C
<b>Ignition temperature:</b>	425	°C.
<b>Danger of explosion:</b>	Lower limit of explosion	3.5 % Kühn-
	Upper limit of explosion	15 % Birett
<b>Density:</b>	ca. 0.8	g/cm <sup>3</sup>
<b>Vapor pressure:</b>	n.d.	
<b>pH:</b>	neutral	
<b>Solubility in/miscibility with:</b>		
<b>Water:</b>	soluble	
<b>Content of solvents:</b>		
<b>Organic solvents:</b>	ethanol	

### 10. Stability and reactivity

**Thermal decomposition / conditions to be avoided**

**Possibility of hazardous reactions:** no information available

**Conditions to avoid:** Heat, direct sunlight and radical starters. The product is stable by proper use. Avoid contact with alkali metals, hydrides, strong oxidizing agents, radical initiators.

**Hazardous decomposition products:** no information available.

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
Printing date: 27.02.2017

Version: 1

### **11. Toxicological information**

**Acute toxicity: (ethanol)**

LD<sub>50</sub> (oral, rat) 6200 mg/kg

**Primary irritation:**

**skin:** degreasing effect

**eye:** irritation

**Sensitization:** no information available

**Additional toxicological information:**

Allergic reactions after long-term exposure are possible.

When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

### **12. Ecological information**

**General information:** Avoid transfer into environment

**Classification of water endangerment:** WGK=1 (slightly polluting substance).

### **13. Disposal considerations**

**Product:**

**Recommendation:**

Use proper landfill disposal or incineration in accordance with local, state and federal regulations.

**Uncleaned packaging:**

**Recommendation:** Disposal must be made according to official regulations

### **14. Transport information**

**Land transport ADR/RID and GGVS/GGVE (Germany):**

UN 1170 ETHANOL, 3, II

**Sea transport IMDG:**

UN 1170 ETHANOL, 3, II

EmS F-E S-D

**Air transport ICAO-TI and IATA-DGR:**

UN 1170 ETHANOL, 3, II

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.

### **15. Regulatory information**

**Classification according to EC-guidelines:**

The product is a medical device according to EC-directive 93/42 EEC.

**Chemical safety assessment:**

A Chemical Safety Assessment has not been carried out.

### **16. Other information:**

**Changes compared with the previous version:**

Adaptation according to 1907/2006/EG, Article 31

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

### **1. Identification of the substance/mixture and of the company/undertaking**

**Product Identifier**

**Trade name:** VITA ADIVA F-CEM

**Use of the preparation:** Dental product

**Details of the supplier of the safety data sheet:**

**Manufacturer/supplier:**

Harvard Dental International GmbH,  
 Margaretenstraße 2-4,  
 15366 Hoppegarten, Germany  
 info@harvard-dental-international.de

**Emergency telephone number:**

+ 49 (0) 30 99 28 978-0 (Quality Management)

### **2. Hazards identification**

**Classification of the substance or mixture:**

**Classification according to Regulation (EC) No 1272/2008:**

Skin Irrit. 2, H315 Causes skin irritation

Eye Irrit. 2, H319 Causes serious eye irritation

Skin Sens. 1; H317 May cause allergic skin reaction

**Information concerning particular hazards for human and environment:**

The product has to be labeled due to calculation procedure for the "General Classification guideline for preparations of the EU" in the latest version.

**Label elements**

**Labeling according to Regulation (EC) 1272/2008:**

**Hazard pictograms**



**Signal word**

Warning

**Hazard statements**

H315 Causes skin irritation

H319 Causes serious eye irritation

H317 May cause allergic skin reaction

**Precautionary statements**

P280 Wear protective gloves/protective clothing/eye protection/face protection

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P302+P352 IF ON SKIN: Wash with plenty of soap and water

**Labeling according to EU guidelines:**

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

**Other hazards**

**Results of PBT and vPvB assessment**

**PBT:** not applicable

**vPvB:** not applicable

## Safety Data Sheet

### according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

#### **3. Composition/Information on Ingredients**

**Chemical characterization:**

**Description:**

Mixture of resin based on Bis-GMA, catalyst, stabilizer, pigments

**Hazardous components:**

Methacrylates                      GHS07, H315, H319, H317, H335

#### **4. First-aid measures**

**General information:** Seek medical advice, if symptoms occur, which can be caused by the product

**after skin contact:** wash with plenty of water and soap

**after inhalation:** fresh air

**after eye contact:** rinse with plenty of water and consult a doctor

**after swallowing:** If symptoms persist consult a doctor

#### **5. Fire fighting measures**

**Extinguishing media:**

CO<sub>2</sub>, powder, water spray.

Use extinguishing method and media suitable to surrounding fire.

**Special hazards arising from the substance or mixture** No further relevant information available

**Protective equipment** Wear self-contained respiratory protective equipment

#### **6. Accidental release measures**

**Personal precautions:** not required

**Environmental precautions:**

Do not allow to enter drains/surface water/ ground-water.

**Methods for cleaning up:**

Pick up mechanically. Clean up affected area.

**Reference to other sections:**

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment

See Section 13 for disposal information

#### **7. Handling and storage**

**Handling:**

This product should only be supplied or handled to dentists or dental technicians or upon their instructions.

**Information for safe handling:** Observe normal care for working with chemicals.

**Information about fire and explosion protection:** No special measures required

**Storage:** Store at a dry place. Store tightly closed. Storage temperature <25°C.

**Requirements at storerooms and containers:** Store only in the original containers

**Further information on storage conditions:**

Store away from food stuffs.

**Specific end uses:** No further relevant information available

## Safety Data Sheet

### according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017

Version: 1

Printing date: 27.02.2017

#### 8. Exposure controls/personal protection

##### Components with limits of values to be supervised at the workplace:

The product contains no relevant quantities of components with limits of values to be supervised at the workplace.

**Additional information:** The lists valid during the making were used as basis

##### Personal protective equipment:

##### General protective and hygienic measures:

The usual precautionary measures are to be adhered to when handling chemicals.

Keep away from foodstuffs and beverages. Avoid contact with eyes or skin

##### Hand protection: Protective gloves.

Due to missing tests no recommendation to the glove material can be given for the product.

**Eye protection:** Safety glasses

**Respiratory protection:** Required, when aerosols are generated

#### 9. Physical and chemical properties

##### Information on basic physical and chemical properties

##### Appearance

<b>Form:</b>	flowable paste
<b>Color:</b>	tooth colored
<b>Odor:</b>	odorless
	<u>Value                  Unit                  Method</u>

##### Safety relevant basic data

**Melting point/melting range:** n.a.

**Boiling point/boiling range:** n.a.

**Flash point:** n.d.

**Ignition temperature:** n.d.

**Danger of explosion:** none

**Density:** 1.72 (20°C) g/cm<sup>3</sup>

**Vapor pressure:** n.a.

**pH:** neutral

**Solubility in/miscibility with:** Partially soluble in toluene, silicone oil

**Water:** insoluble

**Content of solvents:** none

**Organic solvents:** none

#### 10. Stability and reactivity

##### Thermal decomposition / conditions to be avoided

**Possibility of hazardous reactions:** no information available

**Conditions to avoid:** Heat, direct sunlight and radical starters. The product is stable by proper use. Avoid contact with alkali metals, hydrides, strong oxidizing agents, radical initiators.

**Hazardous decomposition products:** No information available.

**Safety Data Sheet**  
**according to 1907/2006/EC, Article 31**

Date of issue: 01.02.2017  
Printing date: 27.02.2017

Version: 1

<b>11. Toxicological information</b> <b>Acute toxicity:</b> no data available <b>Primary irritation:</b> <b>skin:</b> irritant <b>eye:</b> irritant <b>Sensitization:</b> sensitization by skin contact possible <b>Additional toxicological information:</b> Allergic reactions to the product after long-term exposure are possible. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.
<b>12. Ecological information</b> <b>General information:</b> Avoid transfer into environment <b>Classification of water endangerment:</b> WGK=1 (slightly polluting substance).
<b>13. Disposal considerations</b> <b>Product:</b> <b>Recommendation:</b> Small quantities can be disposed of with household waste. Use proper landfill disposal or incineration in accordance with local, state and federal regulations. <b>Uncleaned packaging:</b> <b>Recommendation:</b> Disposal must be made according to official regulations
<b>14. Transport information</b> <b>Land transport ADR/RID and GGVS/GGVE (Germany):</b> no subject to transport regulations <b>Sea transport IMDG:</b> no subject to transport regulations <b>Air transport ICAO-TI and IATA-DGR:</b> no subject to transport regulations
<b>15. Regulatory information</b> <b>Classification according to EC-guidelines:</b> The product is a medical device according to EC-directive 93/42 EEC. <b>Chemical safety assessment:</b> A Chemical Safety Assessment has not been carried out.
<b>16. Other information:</b> <b>Changes compared with the previous version:</b> Adaptation according to 1907/2006/EG, Article 31  The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 1 of 7

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

VITA ADIVA OXY-PREVENT

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

###### Use of the substance/mixture

Auxiliary for manufacture of dental products

##### 1.3. Details of the supplier of the safety data sheet

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
e-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
e-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	

##### 1.4. Emergency telephone number:

+49-(0)761-19240

##### Further Information

medical device

#### SECTION 2: Hazards identification

##### 2.1. Classification of the substance or mixture

###### UN GHS (ST/SG/AC.10/11/Rev.8)

Hazard categories:

Flammable liquid: Flam. Liq. 4

Hazard Statements:

Combustible liquid.

##### 2.2. Label elements

###### UN GHS (ST/SG/AC.10/11/Rev.8)

Signal word: Warning

###### Hazard statements

H227 Combustible liquid.

###### Precautionary statements

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P403 Store in a well-ventilated place.

##### 2.3. Other hazards

No information available.

#### SECTION 3: Composition/information on ingredients

##### 3.2. Mixtures

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 2 of 7

#### Chemical characterization

Substance, organic

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### After inhalation

Provide fresh air.

#### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

#### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

#### After ingestion

Rinse mouth immediately and drink plenty of water.

### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

### 5.2. Special hazards arising from the substance or mixture

Non-flammable. Vapours can form explosive mixtures with air.

### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

#### Additional information

Use water spray jet to protect personnel and to cool endangered containers. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### General advice

Use personal protection equipment.

### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

### 6.3. Methods and material for containment and cleaning up

#### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## SECTION 7: Handling and storage

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 3 of 7

#### 7.1. Precautions for safe handling

##### **Advice on safe handling**

No special measures are necessary.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

##### **Advice on general occupational hygiene**

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat, drink, smoke, sniff.

#### 7.2. Conditions for safe storage, including any incompatibilities

##### **Requirements for storage rooms and vessels**

Keep container tightly closed.

##### **Hints on joint storage**

No special measures are necessary.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

##### **Exposure limits (EH40)**

CAS No	Substance	ppm	mg/m <sup>3</sup>	fibres/ml	Category	Origin
56-81-5	Glycerol, mist	-	10		TWA (8 h)	WEL

#### 8.2. Exposure controls

##### **Individual protection measures, such as personal protective equipment**

##### **Eye/face protection**

Wear eye/face protection.

##### **Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. NBR (Nitrile rubber)

##### **Skin protection**

Use of protective clothing.

##### **Respiratory protection**

In case of inadequate ventilation wear respiratory protection. Open windows to ensure natural ventilation.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state: Liquid  
 Colour:  
 Odour: characteristic

##### **Changes in the physical state**

Melting point/freezing point: not determined  
 Boiling point or initial boiling point and boiling range: 192 °C  
 Flash point: 180 °C

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 4 of 7

#### Flammability

Solid/liquid: not applicable  
Gas: not applicable

#### Explosive properties

The product is not: Explosive.

#### Self-ignition temperature

Solid: not applicable  
Gas: not applicable

Decomposition temperature: not determined

pH-Value: not determined

#### Solubility in other solvents

not determined

Partition coefficient n-octanol/water: not determined

Vapour pressure: not determined

Density: not determined

Relative vapour density: not determined

#### 9.2. Other information

##### Information with regard to physical hazard classes

Oxidizing properties  
Not oxidising.

##### Other safety characteristics

Solid content: not determined

Evaporation rate: not determined

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

#### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

#### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

#### 10.4. Conditions to avoid

none

#### 10.5. Incompatible materials

No information available.

#### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

### SECTION 11: Toxicological information

#### 11.1. Information on hazard classes as defined in GB CLP Regulation

##### Acute toxicity

Based on available data, the classification criteria are not met.

##### Irritation and corrosivity

Based on available data, the classification criteria are not met.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 5 of 7

#### Sensitising effects

Based on available data, the classification criteria are not met.

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

### SECTION 12: Ecological information

#### 12.1. Toxicity

The product is not: Ecotoxic.

#### 12.2. Persistence and degradability

The product has not been tested.

#### 12.3. Bioaccumulative potential

The product has not been tested.

#### 12.4. Mobility in soil

The product has not been tested.

#### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### 12.7. Other adverse effects

No information available.

#### Further information

Avoid release to the environment.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Wash with plenty of water. Completely emptied packages can be recycled.

### SECTION 14: Transport information

#### Marine transport (IMDG)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Air transport (ICAO-TI/IATA-DGR)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 6 of 7

**14.4. Packing group:** No dangerous good in sense of this transport regulation.

**14.5. Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

**14.6. Special precautions for user**

No information available.

**14.7. Maritime transport in bulk according to IMO instruments**

not applicable

## SECTION 15: Regulatory information

### National regulatory information

## SECTION 16: Other information

### Changes

This data sheet contains changes from the previous version in section(s): 9.

### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 7 of 7

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 1 of 7

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

VITA ADIVA OXY-PREVENT

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

###### Use of the substance/mixture

Auxiliary for manufacture of dental products

##### 1.3. Details of the supplier of the safety data sheet

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
e-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
e-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	

##### 1.4. Emergency telephone number:

+49-(0)761-19240

##### Further Information

medical device

#### SECTION 2: Hazards identification

##### 2.1. Classification of the substance or mixture

###### UN GHS (ST/SG/AC.10/11/Rev.8)

Hazard categories:

Flammable liquid: Flam. Liq. 4

Hazard Statements:

Combustible liquid.

##### 2.2. Label elements

###### UN GHS (ST/SG/AC.10/11/Rev.8)

Signal word: Warning

###### Hazard statements

H227 Combustible liquid.

###### Precautionary statements

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P403 Store in a well-ventilated place.

##### 2.3. Other hazards

No information available.

#### SECTION 3: Composition/information on ingredients

##### 3.2. Mixtures

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 2 of 7

#### Chemical characterization

Substance, organic

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### After inhalation

Provide fresh air.

##### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

##### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

##### After ingestion

Rinse mouth immediately and drink plenty of water.

#### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

##### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

#### 5.2. Special hazards arising from the substance or mixture

Non-flammable. Vapours can form explosive mixtures with air.

#### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

#### Additional information

Use water spray jet to protect personnel and to cool endangered containers. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

##### General advice

Use personal protection equipment.

#### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

#### 6.3. Methods and material for containment and cleaning up

##### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

### SECTION 7: Handling and storage

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 3 of 7

#### 7.1. Precautions for safe handling

##### Advice on safe handling

No special measures are necessary.

##### Advice on protection against fire and explosion

No special fire protection measures are necessary.

##### Advice on general occupational hygiene

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat, drink, smoke, sniff.

#### 7.2. Conditions for safe storage, including any incompatibilities

##### Requirements for storage rooms and vessels

Keep container tightly closed.

##### Hints on joint storage

No special measures are necessary.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

##### Exposure limits (EH40)

CAS No	Substance	ppm	mg/m <sup>3</sup>	fibres/ml	Category	Origin
56-81-5	Glycerol, mist	-	10		TWA (8 h)	WEL

#### 8.2. Exposure controls

##### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. NBR (Nitrile rubber)

##### Skin protection

Use of protective clothing.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Open windows to ensure natural ventilation.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state: Liquid  
 Colour:  
 Odour: characteristic

##### Changes in the physical state

Melting point/freezing point: not determined  
 Boiling point or initial boiling point and boiling range: 192 °C  
 Flash point: 180 °C

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 4 of 7

#### Flammability

Solid/liquid: not applicable  
Gas: not applicable

#### Explosive properties

The product is not: Explosive.

#### Self-ignition temperature

Solid: not applicable  
Gas: not applicable

Decomposition temperature: not determined

pH-Value: not determined

#### Solubility in other solvents

not determined

Partition coefficient n-octanol/water: not determined

Vapour pressure: not determined

Density: not determined

Relative vapour density: not determined

#### 9.2. Other information

##### Information with regard to physical hazard classes

Oxidizing properties  
Not oxidising.

##### Other safety characteristics

Solid content: not determined

Evaporation rate: not determined

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

#### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

#### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

#### 10.4. Conditions to avoid

none

#### 10.5. Incompatible materials

No information available.

#### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

### SECTION 11: Toxicological information

#### 11.1. Information on hazard classes as defined in GB CLP Regulation

##### Acute toxicity

Based on available data, the classification criteria are not met.

##### Irritation and corrosivity

Based on available data, the classification criteria are not met.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 5 of 7

#### Sensitising effects

Based on available data, the classification criteria are not met.

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

### SECTION 12: Ecological information

#### 12.1. Toxicity

The product is not: Ecotoxic.

#### 12.2. Persistence and degradability

The product has not been tested.

#### 12.3. Bioaccumulative potential

The product has not been tested.

#### 12.4. Mobility in soil

The product has not been tested.

#### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### 12.7. Other adverse effects

No information available.

#### Further information

Avoid release to the environment.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Wash with plenty of water. Completely emptied packages can be recycled.

### SECTION 14: Transport information

#### Marine transport (IMDG)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Air transport (ICAO-TI/IATA-DGR)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 6 of 7

**14.4. Packing group:** No dangerous good in sense of this transport regulation.

**14.5. Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

**14.6. Special precautions for user**

No information available.

**14.7. Maritime transport in bulk according to IMO instruments**

not applicable

## SECTION 15: Regulatory information

### National regulatory information

## SECTION 16: Other information

### Changes

This data sheet contains changes from the previous version in section(s): 9.

### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA ADIVA OXY-PREVENT

Revision date: 02.06.2022

Product code: 308

Page 7 of 7

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
Printing date: 27.02.2017

Version: 1

### **1. Identification of the substance/mixture and of the company/undertaking**

#### **Product Identifier**

**Trade name:** VITA ADIVA TE-CEM

**Use of the preparation:** Dental product

**Details of the supplier of the safety data sheet:**

#### **Manufacturer/supplier:**

Harvard Dental International GmbH,  
Margaretenstraße 2-4,  
15366 Hoppegarten, Germany  
info@harvard-dental-international.de

#### **Emergency telephone number:**

+ 49 (0) 30 99 28 978-0 (Quality Management)

### **2. Hazards identification**

#### **Classification of the substance or mixture:**

#### **Classification according to Regulation (EC) No 1272/2008:**

Skin Irritation, Category 2; H315 Causes skin irritation

Eye Irritation, Category 2; H319 Causes serious eye irritation

STOT SE, Category 3, H335 May cause respiratory irritation

Acute aquatic toxicity, Category 1, H400 Very toxic to aquatic life

Chronic aquatic toxicity, Category 1, H410 Very toxic to aquatic life with long lasting effects

#### **Information concerning particular hazards for human and environment:**

The product has to be labeled due to the calculation procedure of the „General classification guideline for preparations of the EU“ in the latest valid version.

#### **Label elements:**

#### **Labeling according to Regulation (EC) 1272/2008:**

#### **Hazard pictograms**



#### **Signal word**

Warning

#### **Hazard statements**

H315 Causes skin irritation

H319 Causes serious eye irritation

H335 May cause respiratory irritation

H400 Very toxic to aquatic life

#### **Precautionary statements**

P262 Do not get in eyes, on skin or on clothing.

#### **Labeling according to EU guidelines:**

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

#### **Other hazards**

#### **Results of PBT and vPvB assessment**

**PBT:** not applicable

**vPvB:** not applicable

## Safety Data Sheet

### according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

#### 3. Composition/Information on Ingredients

**Chemical characterization:** Mixture

**Description of the mixture:** Mixture of organic acids, zinc oxide and methacrylates

**Hazardous components:**

Methacrylates [Index-No. 607-134-00-4]                      GHS07, H315, H319, H335

Zinc oxide [CAS-No. 1314-13-2]                              GHS09, H410, H411

#### 4. First-aid measures

**General information:** no special measures required

**after skin contact:** wash with plenty of water and soap

**after eye contact:** Rinse opened eye for several minutes under running water and consult a doctor

**after swallowing:** Seek medical advice immediately

#### 5. Fire fighting measures

**Extinguishing media:**

CO<sub>2</sub>, powder, water spray.

Use extinguishing method and media suitable to surrounding fire

**Special hazards arising from the substance or mixture** No further relevant information available

**Protective equipment** Wear self-contained respiratory protective equipment

#### 6. Accidental release measures

**Personal precautions, protective equipment and emergency procedures**

**Environmental precautions:**

Do not allow to enter drains/surface water/ ground-water.

**Methods for cleaning up:**

Take up mechanically. Clean up affected area. Forward for disposal.

**Reference to other sections:**

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment

See Section 13 for disposal information

#### 7. Handling and storage

**Handling:**

This product should only be supplied or handled to dentists or dental technicians or upon their instructions.

**Precautions for safe handling:** Observe normal care for working with chemicals

**Information about fire and explosion protection:** No special measures required

**Storage:**

**Requirements at storerooms and containers:** Store only in the original containers

**Further information on storage conditions:**

Store at dry place, storage temperature < 25°C. Store away from food stuffs.

**Specific end uses:** No further relevant information available

## Safety Data Sheet

### according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

#### 8. Exposure controls/personal protection

**Components with limits of values to be supervised at the workplace:**

The product contains no relevant quantities of components with limit values to be supervised at the workplace.

**Additional information:** The lists valid during the making were used as basis

**Personal protective equipment:**

**General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.  
 Keep away from foodstuffs and beverages.

**Hand protection**

Protective gloves.

Due to missing tests no recommendation to the glove material can be given for the product.

**Eye protection:** safety glasses

**Respiratory protection:** Required, when aerosols are generated

#### 9. Physical and chemical properties

**Information on basic physical and chemical properties**

**Appearance**

<b>Form:</b>	paste		
<b>Color:</b>	Tooth colored		
<b>Odor:</b>	odorless		
	<u>Value</u>	<u>Unit</u>	<u>Method</u>

**Safety relevant basic data**

<b>Melting point/melting range:</b>	n.a.		
<b>Boiling point/boiling range:</b>	n.a.		
<b>Flash point:</b>	n.a.		
<b>Ignition temperature:</b>	Product is not selfigniting		
<b>Danger of explosion:</b>	None		
<b>Density:</b>	1,56 (20 °C)	g/cm <sup>3</sup>	
<b>Vapor pressure:</b>	Not applicable (non volatile)		
<b>Solubility in/miscibility with:</b>	slightly in toluene, organic solvents		
<b>Water:</b>	insoluble		
<b>Content of solvents:</b>	none		
<b>Organic solvents:</b>	none		

#### 10. Stability and reactivity

**Thermal decomposition / conditions to be avoided**

No decomposition if used according to specifications.

**Possibility of hazardous reactions:** No dangerous reaction known

**Conditions to avoid:** avoid heat, direct sunlight

**Incompatible materials:** Peroxides, metal powders, strong oxidizers, free radical initiators

**Hazardous decomposition products:** No dangerous decomposition products known.

#### 11. Toxicological information

Allergic reactions to methacrylates by chronic exposure to the product are possible

**Acute toxicity:**

**Primary irritation:**

**skin:** slightly irritating

**eye:** irritating effect

**Sensitization:** no sensitizing effect known

## Safety Data Sheet

### according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017

Version: 1

Printing date: 27.02.2017

**Additional toxicological information:**

When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

**12. Ecological information****Toxicity:**

**Aquatic toxicity:** May cause long term adverse effect in the aquatic environment

**Persistence and degradability:** No further relevant information available

**Behavior in environmental systems:**

**Bioaccumulative potential** No further relevant information available

**Mobility in soil** No further relevant information available

**General information:** Do not allow product into surface water or sanitary sewer system

**13. Disposal considerations****Product:****Recommendation:**

Use proper landfill disposal or incineration in accordance with local, state and federal regulations

**Uncleaned packaging:** Handle contaminated packaging in the same way as the substance itself.

**Recommendation:** Disposal must be made according to official regulations

**14. Transport information****Land transport ADR/RID and GGVS/GGVE (Germany)**

no subject to transport regulations.

**Sea transport IMDG**

no subject to transport regulations.

**Air transport ICAO-TI and IATA-DGR**

no subject to transport regulations.

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.

**15. Regulatory information****Classification according to EC-guidelines:**

The product is a medical device according to EC-directive 93/42 EEC.

**Chemical safety assessment:**

A Chemical Safety Assessment has not been carried out.

**16. Other information:****Changes compared with the previous version:**

Adaptation according to 1907/2006/EG, Article 31

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

### **1. Identification of the substance/mixture and of the company/undertaking**

**Product Identifier**

**Trade name:** VITA ADIVA T-BOND I

**Use of the preparation:** Dental product

**Details of the supplier of the safety data sheet:**

**Manufacturer/supplier:**

Harvard Dental International GmbH,

Margaretenstraße 2-4,

15366 Hoppegarten, Germany

info@harvard-dental-international.de

**Emergency telephone number:**

+ 49 (0) 30 99 28 978-0 (Quality Management)

### **2. Hazards identification**

**Classification of the substance or mixture:**

**Classification according to Regulation (EC) No 1272/2008:**

Skin Irrit. 2, H315 Causes skin irritation

Eye Irrit. 2, H319 Causes serious eye irritation

STOT SE 3, H335 May cause respiratory irritation.

**Information concerning particular hazards for human and environment:**

The product has to be labeled due to calculation procedure for the "General Classification guideline for preparations of the EU" in the latest version.

**Label elements**

**Labeling according to Regulation (EC) 1272/2008:**

**Hazard pictograms**



**Signal word**

Danger

**Hazard statements**

H315 Causes skin irritation

H319 Causes serious eye irritation

H335 May cause respiratory irritation

**Precautionary statements**

P280 Wear protective gloves/protective clothing/eye protection/face protection

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P302+P352 IF ON SKIN: Wash with plenty of soap and water

**Labeling according to EU guidelines:**

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

**Other hazards**

**Results of PBT and vPvB assessment**

**PBT:** not applicable

**vPvB:** not applicable

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

### 3. Composition/Information on Ingredients

**Chemical characterization:**

**Description:**

Modified acrylic acid, poly acrylic acid, methacrylates, catalysts, stabilizers in ethanol

**Hazardous components:**

607-134-00-4 (Index No.)	Methacrylates	GHS07, H315, H319, H335
64-17-5 (CAS-No.)	Ethanol	GHS02, H225

### 4. First-aid measures

**General information:** Seek medical advice, if symptoms occur, which can be caused by the product

**after skin contact:** wash with plenty of water and soap

**after inhalation:** fresh air

**after eye contact:** Rinse opened eye for several minutes under running water and consult a doctor

**after swallowing:** Consult a physician.

### 5. Fire fighting measures

**Extinguishing media:**

CO<sub>2</sub>, powder, water spray.

Use extinguishing method and media suitable to surrounding fire.

**Special hazards arising from the substance or mixture** No further relevant information available

**Protective equipment** Wear self-contained respiratory protective equipment

### 6. Accidental release measures

**Personal precautions, protective equipment and emergency procedures:**

Do not inhale vapours/aerosols. Avoid substance contact. Ensure supply of fresh air in enclosed rooms

**Environmental precautions:**

Do not allow to enter drains/surface water/ ground-water.

**Methods for cleaning up:**

Absorb with liquid binding material. Clean up affected area.

**Reference to other sections:**

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment

See Section 13 for disposal information

### 7. Handling and storage

**Handling:**

This product should only be supplied or handled to dentists or dental technicians or upon their instructions.

**Precautions for safe handling:** Observe normal care for working with chemicals

**Information about fire and explosion protection:** Keep away from heat and fire. Take measure to prevent electrostatic charge

**Storage:** Store at a dry place. Store tightly closed. Storage temperature <25°C.

**Requirements at storerooms and containers:** Store only in the original containers

**Further information on storage conditions:**

Store away from food stuffs.

**Specific end uses:** No further relevant information available

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

### 8. Exposure controls/personal protection

**Components with limits of values to be supervised at the workplace:**

CAS-No.	name	art	value	unit
64-17-5	Ethanol	TRGS 900	500	ml/m <sup>3</sup> (ppm)
			or 960	mg/m <sup>3</sup>

**Additional information:** The lists valid during the making were used as basis

**Personal protective equipment:**

**General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

Keep away from foodstuffs and beverages. Avoid contact with eyes or skin

**Hand protection**

Protective gloves.

Due to missing tests no recommendation to the glove material can be given for the product.

**Eye protection:** safety glasses

**Respiratory protection:** Required, when aerosols are generated

### 9. Physical and chemical properties

**Information on basic physical and chemical properties**

**Appearance**

<b>Form:</b>	liquid						
<b>Color:</b>	colorless						
<b>Odor:</b>	characteristic						
	<table border="1"> <thead> <tr> <th>Value</th> <th>Unit</th> <th>Method</th> </tr> </thead> <tbody> <tr> <td> </td> <td> </td> <td> </td> </tr> </tbody> </table>	Value	Unit	Method			
Value	Unit	Method					

**Safety relevant basic data**

**Melting point/melting range:** n.a.

**Boiling point/boiling range:** n.a.

**Flash point:** n.a.

**Ignition temperature:** no information available

**Danger of explosion:**

**Density:** ca. 1.0 g/cm<sup>3</sup>

**Vapor pressure:** n.d.

**pH:** ca. 3

**Solubility in/miscibility with:**

**Water:** soluble

**Content of solvents:**

**Organic solvents:** ethanol

### 10. Stability and reactivity

**Thermal decomposition / conditions to be avoided**

**Possibility of hazardous reactions:** no information available

**Conditions to avoid:** Heat, direct sunlight and radical starters. The product is stable by proper use. Avoid contact with alkali metals, hydrides, strong oxidizing agents, radical initiators.

**Hazardous decomposition products:** No information available.

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
Printing date: 27.02.2017

Version: 1

### **11. Toxicological information**

**Acute toxicity: (ethanol)**

LD<sub>50</sub> (oral, rat) 6200 mg/kg

**Primary irritation:**

**skin:** irritant

**eye:** irritant

**Sensitization:** no information available

**Additional toxicological information:**

Allergic reactions to methacrylates after long-term exposure are possible. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

### **12. Ecological information**

**General information:** Avoid transfer into environment

**Classification of water endangerment:** WGK=1 (slightly polluting substance).

### **13. Disposal considerations**

**Product:**

**Recommendation:**

Use proper landfill disposal or incineration in accordance with local, state and federal regulations.

**Uncleaned packaging:**

**Recommendation:** Disposal must be made according to official regulations

### **14. Transport information**

**Land transport ADR/RID and GGVS/GGVE (Germany):**

No subject to transport regulations

**Sea transport IMDG:**

No subject to transport regulations

**Air transport ICAO-TI and IATA-DGR:**

No subject to transport regulations

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.

### **15. Regulatory information**

**Classification according to EC-guidelines:**

The product is a medical device according to EC-directive 93/42 EEC.

**Chemical safety assessment:**

A Chemical Safety Assessment has not been carried out.

### **16. Other information:**

**Changes compared with the previous version:**

Adaptation according to 1907/2006/EG, Article 31

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
Printing date: 27.02.2017

Version: 1

### 1. Identification of the substance/mixture and of the company/undertaking

#### Product Identifier

**Trade name: VITA ADIVA ZR-PRIME**

**Use of the preparation:** Dental product

**Details of the supplier of the safety data sheet:**

#### Manufacturer/supplier:

Harvard Dental International GmbH,  
Margaretenstraße 2-4,  
15366 Hoppegarten, Germany  
info@harvard-dental-international.de

#### Emergency telephone number:

+ 49 (0) 30 99 28 978-0 (Quality Management)

### 2. Hazards identification

#### Classification of the substance or mixture:

#### Classification according to Regulation (EC) No 1272/2008:

Flam. Liq. 2, H225 Highly flammable liquid and vapour

Eye Irrit. 2, H319 Causes serious eye irritation

STOT SE 3, H336 May cause drowsiness or dizziness.

#### Information concerning particular hazards for human and environment:

The product has to be labeled due to calculation procedure for the "General Classification guideline for preparations of the EU" in the latest version.

#### Label elements

#### Labeling according to Regulation (EC) 1272/2008:

#### Hazard pictograms



#### Signal word

Danger

#### Hazard statements

H225 Highly flammable liquid and vapour

H319 Causes serious eye irritation

H336 May cause drowsiness or dizziness

#### Precautionary statements

P210 Keep away from heat/sparks/open flames/hot surfaces – No smoking.

P233 Keep container tightly closed..

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

#### Labeling according to EU guidelines:

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

#### Other hazards

#### Results of PBT and vPvB assessment

**PBT:** not applicable

**vPvB:** not applicable

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

<b>3. <u>Composition/Information on Ingredients</u></b>
<p><b>Chemical characterization:</b>  <b>Beschreibung:</b>          Preparation based on acetone</p> <p><b>Hazardous components:</b>          CAS 67-64-1(EINECS: 200-662-2) Acetone GHS02, GHS07, H225, H319, H336</p>
<b>4. <u>First-aid measures</u></b>
<p><b>General information:</b> Seek medical advice, if symptoms occur, which can be caused by the product</p> <p><b>after skin contact:</b> Wash with plenty of water and soap</p> <p><b>after inhalation:</b> Fresh air. Consult a physician if feeling unwell.</p> <p><b>after eye contact:</b> Rinse opened eye for several minutes under running water and consult an ophthalmologist</p> <p><b>after swallowing:</b> Caution: Aspiration hazard. Keep respiratory tract clear. Call a physician.</p>
<b>5. <u>Fire fighting measures</u></b>
<p><b>Extinguishing media:</b>          CO<sub>2</sub>, powder, water spray.          Use extinguishing method and media suitable to surrounding fire.</p> <p><b>Special hazards arising from the substance or mixture</b> No further relevant information available</p> <p><b>Protective equipment</b> Wear self-contained respiratory protective equipment</p>
<b>6. <u>Accidental release measures</u></b>
<p><b>Personal precautions, protective equipment and emergency procedures:</b>          Do not inhale vapours/aerosols. Avoid substance contact. Ensure supply of fresh air in enclosed rooms</p> <p><b>Environmental precautions:</b>          Do not allow to enter drains/surface water/ ground-water.</p> <p><b>Methods for cleaning up:</b>          Absorb with liquid binding material. Clean up affected area. Ensure supply of fresh air.</p> <p><b>Reference to other sections:</b>          See Section 7 for information on safe handling          See Section 8 for information on personal protection equipment          See Section 13 for disposal information</p>
<b>7. <u>Handling and storage</u></b>
<p><b>Handling:</b>          This product should only be supplied or handled to dentists or dental technicians or upon their instructions.</p> <p><b>Precautions for safe handling:</b> Observe normal care for working with chemicals</p> <p><b>Information about fire and explosion protection:</b> Keep away from heat and fire. Take measure to prevent electrostatic charge</p> <p><b>Storage:</b> Store at a dry place. Store tightly closed. Storage temperature &lt;25°C.</p> <p><b>Requirements at storerooms and containers:</b> Store only in the original containers</p> <p><b>Further information on storage conditions:</b>          Store away from food stuffs.</p> <p><b>Specific end uses:</b> No further relevant information available</p>

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017

Version: 1

Printing date: 27.02.2017

### 8. Exposure controls/personal protection

**Components with limits of values to be supervised at the workplace:**

CAS-No.	name	art	value	unit
67-64-1	Acetone	TRGS 900	500	ml/m <sup>3</sup> (ppm)
			or 1200	mg/m <sup>3</sup>

**Additional information:** The lists valid during the making were used as basis

**Personal protective equipment:**

**General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

Keep away from foodstuffs and beverages. Avoid contact with eyes or skin

**Hand protection**

Protective gloves.

Due to missing tests no recommendation to the glove material can be given for the product.

**Eye protection:** safety glasses

**Respiratory protection:** Required, when aerosols are generated

### 9. Physical and chemical properties

**Information on basic physical and chemical properties**

**Appearance**

<b>Form:</b>	liquid						
<b>Color:</b>	colorless						
<b>Odor:</b>	characteristic						
	<table border="1"> <thead> <tr> <th>Value</th> <th>Unit</th> <th>Method</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> <td></td> </tr> </tbody> </table>	Value	Unit	Method			
Value	Unit	Method					

**Safety relevant basic data**

**Melting point/melting range:** n.d.

**Boiling point/boiling range:** 55 °C

**Flash point:** -20 °C DIN 51794

**Ignition temperature:** no information available

**Danger of explosion:**

**Density:** ca. 0.8 g/cm<sup>3</sup>

**Vapor pressure:** n.d.

**Solubility in/miscibility with:**

**Water:** soluble

**Content of solvents:**

**Organic solvents:** Acetone

### 10. Stability and reactivity

**Thermal decomposition / conditions to be avoided**

**Possibility of hazardous reactions:** no information available

**Conditions to avoid:** Heat, sunlight, peroxides, alcali hydroxides, halogens, halogenated hydrocarbons, strong oxidizers, free radical initiators.

**Hazardous decomposition products:** No information available.

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 01.02.2017  
 Printing date: 27.02.2017

Version: 1

### **11. Toxicological information**

**Acute toxicity: (acetone)**

LD<sub>50</sub> (oral, rat) 5800 mg/kg

**Primary irritation:**

**skin:** irritant. Degreasing effect on the skin, possibly followed by secondary inflammation.

**eye:** irritant

**Sensitization:** no information available

**Additional toxicological information:**

Allergic reactions after long-term exposure are possible. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

### **12. Ecological information**

**General information:** Avoid transfer into environment

**Classification of water endangerment:** WGK=1 (slightly polluting substance).

### **13. Disposal considerations**

**Product:**

**Recommendation:**

Use proper landfill disposal or incineration in accordance with local, state and federal regulations.

**Uncleaned packaging:**

**Recommendation:** Disposal must be made according to official regulations

### **14. Transport information**

**Land transport ADR/RID and GGVS/GGVE (Germany):**

UN 1090 ACETON, 3, II

**Sea transport IMDG:**

UN 1090 ACETONE, 3 II

EmS: F-E S-D

**Air transport ICAO-TI and IATA-DGR:**

UN 1090 ACETONE, 3, II

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.

### **15. Regulatory information**

**Classification according to EC-guidelines:**

The product is a medical device according to EC-directive 93/42 EEC.

**Chemical safety assessment:**

A Chemical Safety Assessment has not been carried out.

### **16. Other information:**

**Changes compared with the previous version:**

Adaptation according to 1907/2006/EG, Article 31

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.