

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name

BegoForm (REF 52785)

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses of the substance or mixture

Manufacturing of dental prosthesis in a dental laboratory

Uses advised against

No data available.

1.3 Details of the supplier of the safety data sheet

Address

BEGO Bremer Goldschlägerei

Wilh. Herbst GmbH & Co. KG

Wilhelm-Herbst-Str. 1

28359 Bremen

Telephone no. +49/ 421/ 2028 – 0

Fax no. +49/ 421/ 2028 – 115

e-mail msds@bego.com

Information provided by / telephone

Research & Development Department - Materials, alloys and ceramics; +49/ 421/ 2028 – 130 (Chief Development Officer alloys)

Advice on Safety Data Sheet

msds@bego.com

1.4 Emergency telephone number

For medical advice (in German and English):

+49 (0)551 192 40 (Giftinformationszentrum Nord)

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 (CLP)

STOT RE 2; H373i

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008 (CLP Regulation)

Hazard pictograms



GHS08

Signal word

Warning

Hazardous component(s) to be indicated on label:

Quartz (SiO₂)

Hazard statements

H373i

May cause damage to organs through prolonged or repeated exposure if inhaled.

Precautionary statements

P260

Do not breathe dust.

P314

Get medical advice/attention if you feel unwell.

P403+P233

Store in a well-ventilated place. Keep container tightly closed.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

2.3 Other hazardsPBT assessment
No data available.vPvB assessment
No data available.**SECTION 3: Composition/information on ingredients****3.1 Substances**

Not applicable. The product is not a substance.

3.2 Mixtures**Hazardous ingredients**

No	Substance name	Additional information	
	CAS / EC / Index / REACH no	Classification (EC) 1272/2008 (CLP)	Concentration %
1	Quartz (SiO₂)		
	14808-60-7 238-878-4 -	STOT RE 2; H373i	> 10.00 - < 50.00 %b.w.
2	cristobalite		
	14464-46-1 238-455-4 -	STOT RE 2; H373i	> 5.00 - < 20.00 %b.w.
3	zinc oxide		
	1314-13-2 215-222-5 030-013-00-7 01-2119463881-32	Aquatic Acute 1; H400 Aquatic Chronic 1; H410	< 0.50 %b.w.
4	aluminium oxide		
	1344-28-1 215-691-6 -	-	> 10.00 - < 30.00 %b.w.
5	MAGNESIUM OXIDE		
	1309-48-4 215-171-9 -	-	> 10.00 - < 30.00 %b.w.

Full Text for all H-phrases and EUH-phrases: pls. see section 16

No	Note	Specific concentration limits	M-factor (acute)	M-factor (chronic)
3	-	-	M = 1	-

No	Route, target organ, concrete effect
1	H373i inhalational; -; -
2	H373i inhalational; -; -

SECTION 4: First aid measures**4.1 Description of first aid measures****General information**

If the patient is likely to become unconscious, place and transport in stable sideways position. In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air. Irregular breathing/no breathing: artificial respiration. Call a doctor immediately.

After skin contact

In case of contact with skin wash off with water.

After eye contact

Remove contact lenses, irrigate copiously with clean, fresh water for at least 15 minutes, holding the eyelids apart and seek medical advice.

After ingestion

Rinse mouth thoroughly with water. Never give anything by mouth to an unconscious person. Call a doctor immediately.

4.2 Most important symptoms and effects, both acute and delayed

No data available.

4.3 Indication of any immediate medical attention and special treatment needed

No data available.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media

Carbon dioxide; Extinguishing powder; Water spray jet; Foam

Unsuitable extinguishing media

High power water jet

5.2 Special hazards arising from the substance or mixture

No data available.

5.3 Advice for firefighters

Product itself does not burn. Adapt extinguisher and fire-fighting measures to fire in the environment. Use self-contained breathing apparatus. Wear protective clothing.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

Avoid dust formation. Refer to protective measures listed in sections 7 and 8.

For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

6.3 Methods and material for containment and cleaning up

Take up mechanically. Avoid raising dust. When picked up, treat material as prescribed under heading "Disposal considerations".

6.4 Reference to other sections

No data available.

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Keep container tightly closed.

General protective and hygiene measures

Wash hands before breaks and after work. Do not eat, drink or smoke during work time. Remove soiled or soaked clothing immediately. Do not inhale dust. Keep away from foodstuffs and beverages.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

Advice on protection against fire and explosion

No special measures necessary.

7.2 Conditions for safe storage, including any incompatibilities

Technical measures and storage conditions

Keep container tightly closed and dry in a cool, well-ventilated place.

Requirements for storage rooms and vessels

Keep only in the original container.

Advice on storage assembly

Protect against humidity; Do not store together with: explosive substances

7.3 Specific end use(s)

No data available.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

SECTION 8: Exposure controls/personal protection**8.1 Control parameters****Occupational exposure limit values**

No	Substance name	CAS no.	EC no.
1	Dust		
	List of approved workplace exposure limits (WELs) / EH40		
	total inhalable dust		
	TWA	10	mg/m ³
	List of approved workplace exposure limits (WELs) / EH40		
	respirable dust		
	TWA	4	mg/m ³
2	aluminium oxide	1344-28-1	215-691-6
	List of approved workplace exposure limits (WELs) / EH40		
	Aluminium oxides		
	total inhalable dust		
	TWA	10	mg/m ³
	List of approved workplace exposure limits (WELs) / EH40		
	Aluminium oxides		
	respirable dust		
	TWA	4	mg/m ³
3	MAGNESIUM OXIDE	1309-48-4	215-171-9
	List of approved workplace exposure limits (WELs) / EH40		
	Magnesium oxide		
	inhalable dust		
	Mg		
	TWA	10	mg/m ³
	List of approved workplace exposure limits (WELs) / EH40		
	Magnesium oxide		
	fume and respirable dust		
	Mg		
	TWA	4	mg/m ³
4	ammonia, anhydrous	7664-41-7	231-635-3
	List of approved workplace exposure limits (WELs) / EH40		
	Ammonia		
	STEL	25	mg/m ³ 35 ml/m ³
	TWA	18	mg/m ³ 25 ml/m ³
	2000/39/EWG		
	Ammonia, anhydrous		
	STEL	36	mg/m ³ 50 ml/m ³
	TWA	14	mg/m ³ 20 ml/m ³
5	zinc oxide	1314-13-2	215-222-5
	List of approved workplace exposure limits (WELs) / EH40		
	Zinc oxide, fume		
	STEL	10	mg/m ³
	TWA	5	mg/m ³

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

DNEL and PNEC values**DNEL values (worker)**

No	Substance name			CAS / EC no
	Route of exposure	Exposure time	Effect	Value
1	zinc oxide			1314-13-2 215-222-5
	dermal	Long term (chronic)	systemic	83 mg/kg/day
	with reference to: Zn Remarks: insoluble			
	inhalative	Long term (chronic)	systemic	5 mg/m ³
	with reference to: Zn Remarks: insoluble			

DNEL value (consumer)

No	Substance name			CAS / EC no
	Route of exposure	Exposure time	Effect	Value
1	zinc oxide			1314-13-2 215-222-5
	oral	Long term (chronic)	systemic	0.83 mg/kg/day
	with reference to: Zn Remarks: insoluble			
	dermal	Long term (chronic)	systemic	83 mg/kg/day
	with reference to: Zn Remarks: insoluble			
	inhalative	Long term (chronic)	systemic	2.5 mg/m ³
	with reference to: Zn Remarks: insoluble			

PNEC values

No	Substance name		CAS / EC no
	ecological compartment	Type	Value
1	zinc oxide		1314-13-2 215-222-5
	water	fresh water	20.6 µg/l
	with reference to: Zn		
	water	marine water	6.1 µg/l
	with reference to: Zn		
	water	fresh water sediment	117.8 mg/kg
	water	marine water sediment	56.5 mg/kg
	with reference to: Zn, dry weight		
	soil	-	35.6 mg/kg
	with reference to: Zn, dry weight		
	sewage treatment plant	-	100 µg/l

8.2 Exposure controls**Appropriate engineering controls**

No data available.

Personal protective equipment**Respiratory protection**

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of dust formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

Eye / face protection

Safety glasses with side protection shield (EN 166)

Trade name: BegoForm (REF 52785)**Current version :** 2.0.3, issued: 19.10.2015**Replaced version:** 2.0.2, issued: 16.12.2014**Region:** GB**Hand protection**

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product . Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

Other

Normal chemical work clothing.

Environmental exposure controls

No data available.

SECTION 9: Physical and chemical properties**9.1 Information on basic physical and chemical properties**

Form/Colour			
Powder			
white			
Odour			
odourless			
Odour threshold			
No data available			
pH value			
Value	5	-	7
Concentration		80	% (water)

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

Boiling point / boiling range

No data available

Melting point / melting range

No data available

Decomposition point / decomposition range

No data available

Flash point

No data available

Auto-ignition temperature

No data available

Oxidising properties

No data available

Explosive properties

No data available

Flammability (solid, gas)

No data available

Lower flammability or explosive limits

No data available

Upper flammability or explosive limits

No data available

Vapour pressure

No data available

Vapour density

No data available

Evaporation rate

No data available

Relative density

No data available

Density

No data available

Bulk density

Value	1.2	-	1.5	g/cm ³
Reference temperature			21	°C

Solubility in water

No data available

Solubility(ies)

No data available

Partition coefficient: n-octanol/water

No data available

Viscosity

No data available

9.2 Other information**Other information**

The product is hygroscopic.

Product hardens with the addition of water.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

SECTION 10: Stability and reactivity**10.1 Reactivity**

Dangerous reactions are not expected if the product is handled according to its intended use.

10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

10.3 Possibility of hazardous reactions

Dangerous reactions are not to be expected when handling product according to its intended use.

10.4 Conditions to avoid

Protect from humid air and water.

10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products

Under normal heating of the hardened ground-mass in the pre-heat furnace ammonia is released with a concentration of 1 percent by weight. Provide ventilation to exhaust vapours and conduct them directly to the outside.

SECTION 11: Toxicological information**11.1 Information on toxicological effects****Acute oral toxicity**

No data available

Acute dermal toxicity

No data available

Acute inhalational toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Reproduction toxicity

No data available

Carcinogenicity

No data available

STOT-single exposure

No data available

STOT-repeated exposure

No data available

Aspiration hazard

No data available

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Contact with the skin and eyes may cause mechanical irritation. Inhalation of dusts may irritate the respiratory tract. Can cause silicosis.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

SECTION 12: Ecological information**12.1 Toxicity**

Toxicity to fish (acute)
No data available
Toxicity to fish (chronic)
No data available
Toxicity to Daphnia (acute)
No data available
Toxicity to Daphnia (chronic)
No data available
Toxicity to algae (acute)
No data available
Toxicity to algae (chronic)
No data available
Bacteria toxicity
No data available

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

Results of PBT and vPvB assessment	
PBT assessment	No data available.
vPvB assessment	No data available.

12.6 Other adverse effects

No data available.

12.7 Other information

Other information
Do not discharge product unmonitored into the environment.

SECTION 13: Disposal considerations**13.1 Waste treatment methods****Product**

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

Packaging

Residuals must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

SECTION 14: Transport information**14.1 Transport ADR/RID/ADN**

The product is not subject to ADR/RID/ADN regulations.

14.2 Transport IMDG

The product is not subject to IMDG regulations.

14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

14.4 Other information

No data available.

14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

14.6 Special precautions for user

No data available.

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not relevant

SECTION 15: Regulatory information**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulations****Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)**

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

REACH candidate list of substances of very high concern (SVHC) for authorisation

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, PREPARATIONS AND ARTICLES

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances

This product is not subject to Part 1 or 2 of Annex I.

Other regulations

Adhere to the national sanitary and occupational safety regulations when using this product.

15.2 Chemical safety assessment

No data available.

SECTION 16: Other information**Sources of key data used to compile the data sheet:**

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

EC Directives 2000/39/EC, 2006/15/EC, 2009/161/EU

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding chapter.

Trade name: BegoForm (REF 52785)

Current version : 2.0.3, issued: 19.10.2015

Replaced version: 2.0.2, issued: 16.12.2014

Region: GB

Full text of the H- and EUH- phrases drawn up in sections 2 and 3 (provided not already drawn up in these sections)

H400

Very toxic to aquatic life.

H410

Very toxic to aquatic life with long lasting effects.

Department issuing safety data sheet

UMCO Umwelt Consult GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 183, Tel.: +49(40)79 02 36 300, Fax: +49(40)79 02 36 357, e-mail: umco@umco.de

This information is based on our present state of knowledge and experience.

The security data sheet describes products with a view to the security requirements.

However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO Umwelt Consult GmbH.