

Angebot #254905

	Artikel	Stückpreis
 <p>RECIPROC Gutta-Percha ISO 050 • Precise fit to canals shaped with RECIPROC® NiTi instruments • All cold and warm obturation techniques • Phase with low melting point • Injection moulded</p> <p>Endo Easy Efficient®</p>	<p>Packung 60 Stück ISO 050</p> <p>Artikelnummer: 22201</p> <p>ISO Größe <input type="text" value="050"/></p> <p>Verpackungseinheit <input type="text" value="Packung"/></p> <p>Verpackungsinhalt <input type="text" value="60 Stifte"/></p>	<p>20,21 €*[*]</p>